

Vitrile Gloves - Disposable

VARIOUS SIZES AND COLOURS

Used in hospitals, housework, painting and coating industry as well as the printing and dyeing industry and food processing.



Variations, Product Codes & Item Dimensions:

Colour	Product Info	Order Code	Dimensions	Weight
Black	Large	GLOV.V.BL.100.L.I	92 x 38 mm	5.5 g
	Extra Large	GLOV.V.BL.100.XL.I	95 x 37 mm	6 g
Blue	Large	GLOV.V.B.100.L.I	92 x 38 mm	5.5 g
	Extra Large	GLOV.V.B.100.XL.I	95 x 37 mm	6 g

Features:

- Compounded with PVC paste and Nitrile latex
- Softer than vinyl gloves
- Easy to put on when hands are wet
- Odourless
- Non sterile
- Single Use
- Packing: 100pcs/box
- Storage: Store in a cool, dry place (< 40°C) and away from direct sunlight

Be Safe Paramedical cc South Africa

Cri+i Care
 "Progress Saves Lives..."

An ISO 9001:2015 & 13485:2016 Certified Company
 www.be-safe.co.za
 info@be-safe.co.za

CAPE TOWN: +27 21 788 4681
 JOHANNESBURG: +27 11 794 7476
 DURBAN: +27 31 702 3219

Vitrile Gloves - Disposable

VARIOUS SIZES AND COLOURS



Method of use:

- Ensure the gloves fit securely around your fingers and wrist so the gloves don't slip off
- When taking the gloves off, use the inside out method; pinch the cuff at the wrist and pull towards your fingers until the glove comes off. This prevents substances from touching your skin.

Note: Long/sharp nails may damage the gloves

Standards Compliance:

- [ISO 9001:2015](#)
- [ISO 13485:2016](#)
- [FDA](#)
- [CE](#)

Be Safe Paramedical cc South Africa



An ISO 9001:2015 & 13485:2016 Certified Company
www.be-safe.co.za
info@be-safe.co.za

CAPE TOWN: +27 21 788 4681
JOHANNESBURG: +27 11 794 7476
DURBAN: +27 31 702 3219

V01-2020