

The Responder

"Progress Saves Lives..."

Edition 72 | October 2024

From the Editor's Desk

Greetings friends, clients, and fellow healthcare professionals, welcome to the October edition of The Responder.

My goodness Winter really seems to be holding on for dear life. Just when we thought the couple of warm days we've recently had was the real beginnings of Spring, and just like that, we were back to rain and really icy conditions. I guess that's SA for you, four seasons in one day. Having said that, the wet weather certainly could not put a damper on Saturday's Bok game. Well done Bokke on a very convincing win. You once again showed why you are the best team in the world.

We have some exciting news to share with you. We will be exhibiting at the Africa Health Expo, which runs from 22nd to the 24th of October 2024 at the Cape Town International Convention Centre. This is the first time Africa Health is being held in Cape Town and we are really excited to be part of it. We will post our stand numbers on our Social Media pages as soon as we have received them. If you are visiting Africa Health this year, you will find us on the MDMSA and Wesgro pavilions, both of who were instrumental in us being able to exhibit this year. We would love for you to visit our stands to see first-hand, some of the exciting medical devices we manufacture and of course, some of our newest additions. We hope to see you there.



Something I find truly exciting is when we see the products we manufacture being used by Emergency Medical Services in SA. Our CritiCare® brand has become synonymous with quality and durability and it is something we have always prided ourselves in. I would love to share this photo, snapped by Sioux, one of the staff from our Head Office, at O.R. Tambo International Airport recently. To the service, thank you for using and trusting in our products and to the two Medics, thank you for your service.

In this edition of The Responder, our article will focus on the use of tourniquets to control severe external bleeding. As the manufacturer of the popular Rotation Compression Tourniquet (R-C-T), we can say that it is interesting to note that many services are still not using them as often as what one would think they should. Tourniquets, purpose-designed for the control of severe external bleeding on a limb can be truly lifesaving. We would love to see a drive whereby services place tourniquets on all their emergency vehicles and in all their jump bags. I truly believe that it would make a very significant difference to patient care. Something to really think about.

Also, a good time to announce that we will be closing for our year-end holidays on the 13th of December 2024 and reopening on the 6th of January 2025. We strongly suggest you get your orders in to cover you for our break period as early as possible in December. It is a very busy time for us and it is important to us to ensure you receive the stock items you need to cover you for the festive season.

As always, we have some great specials lined up for you, so you will definitely want to check them out for some great savings.

With that said, we hope you enjoy this edition of The Responder.

Yours in Healthcare

Gavin Sutton

The Responder

“Progress Saves Lives...”

Edition 72 | October 2024

Healthcare Article

To Tourniquet or not to Tourniquet... that is the question?

Is it not amazing how certain treatments fall in favour, then out and then back in again? I imagine that when most of us were trained you were taught to err on the side of omission when it came to applying a tourniquet. Effectively, patients would need to be banging on the doors at the Pearly Gates before someone would say, “Hmmm maybe we should try a tourniquet...”

So why were we so jumpy when it came to the application of tourniquets in the good 'ol days? Well from personal experience in my early days as an emergency care provider (which was in the late 80's), I remember my instructor telling me. “Your patient will lose his/her limb”. I remember asking, “Wouldn't it be better to lose a limb than almost your entire blood volume and then your life?” The answer seemed to be in favour of the latter. Now of course this is no reflection on the educator of the day, because the teachings were based on the information that was available at that time. We had no major wars to gather evidence from (which I suppose is good and bad). When peace started going pear-shaped in the Middle East, with mass casualties and some nasty limb injuries (studies began to emerge showing that around 44% of patients who died in the field were as a result of extremity trauma), military medical units had to develop a different set of rules to save lives. And so, the age of the tourniquet reappeared with great success. So how different is the combat environment to a regular Friday night in Cape Town or Johannesburg EMS? Sadly, not very different at all and with that goes a different set of rules.



The RCT

Strangely, many of the medics I come into contact with on a daily basis are still favouring the approach of “as a last resort.” The million-dollar question is, “When does the last resort period start?” Is it just prior to the patient acquiring their harp to begin life outside this world? Or is it when we decide that the patient has a potentially uncontrollable arterial bleed and would benefit from the application of a tourniquet sooner rather than later? I have seen them work, so I am a firm supporter in their use. I also believe that EMS services should be playing a more active role in equipping their crews with this life-saving device. The Philadelphia Police Department recently issued over 5000

tactical tourniquets to their police officers and already there are numerous reports of patients surviving injuries that would no doubt have resulted in the loss of their lives before the tourniquet initiative was implemented. So, the police department is taking the lead where EMS should be setting the standard...



Not for Haemorrhage control

Most tourniquets been used in EMS in South Africa are makeshift and mainly consisted of some form of bandage being twisted around the extremity and held in place by the provider's hands. Trust me; this is no easy



The Responder

"Progress Saves Lives..."

Edition 72 | October 2024

feat at all. In fact, it is really difficult to get any form of pressure to stop flow from the artery. In addition to this, the bandages that were most often used were rather narrow in width and as a result a very narrow margin of pressure was applied to the extremity. We know now that with tourniquets, it is effectively only the medial section of the band from which the majority of pressure is applied.

A true tourniquet (to stop arterial haemorrhaging) should not be confused with the stretchy nylon types we use to site peripheral IV's. They are simply not designed to stem large volume bleeds. They really are only suitable for engorging peripheral veins to site an IV cannula.

Which is the best option and how do we apply it?

Elimination of distal artery flow is the fundamental purpose of a tourniquet. Any tourniquet that cannot stop arterial blood-flow is simply ineffective. One precursor to the tourniquet was the use of tight bandages placed proximal to the point of a traumatic amputation. Most commercially manufactured tourniquets are designed for quick application, usually within 60 seconds. A tourniquet replaces the need for direct pressure, additional bandages or a second care provider during transport. An appropriately trained emergency care provider should be able to effectively stop extremity haemorrhage, accomplish multiple interventions and initiate immediate transport without the need for additional emergency personal

Currently, most medics are trained to place the tourniquet just above the sight of the injury while avoiding placement over a joint. Placing the tourniquet more proximally to the injury on the thickest portion of the extremity helps to limit damage to the underlying tissue and prevent any minor bleeds near the injury. Once positioned, tourniquet tensioning requirements will increase with the size of the limb. The pressure required to halt arterial blood flow is inversely proportional to the width of the tourniquet. Complete occlusion of arterial flow on the lower extremities is extremely difficult, if not impossible, with a 25mm wide tourniquet.

A wider tourniquet will be more effective on a lower extremity, as there is less pressure needed on the greater surface area to successfully occlude arterial blood flow, and less tissue damage will occur.



hand application as well as a Velcro securing strap for extra anchoring.

Once the required effect is achieved, the windlass of the tourniquet should be placed in a docking system so as to avoid it unravelling and releasing the tension initially achieved.



Most purpose designed tourniquets have a width of around 50mm and make use of a windlass or ratchet type tensioner to achieve the effects of arterial occlusion.

Many, such as the Rotation Compression Tourniquet also have a buckle for single hand and dual





The Responder

“Progress Saves Lives...”

Edition 72 | October 2024

How tight should it be applied? Well, the best guide is to look at the haemorrhage site while tensioning up the tourniquet. When haemorrhaging stops, you have achieved enough tension.

How long should it be applied for?

Emergency care providers are cautious when it comes to the use of tourniquets because prolonged use can lead to neurovascular damage and tissue death. We know that tissue death from impaired circulation can occur in as little as two hours. We also know from studies that tourniquets have been left on for over 16 hours without any notable harm.

In the emergency medical service setting, there is a huge focus on response times, limiting scene time and relatively short transfer times to a medical facility. In these instances, it is relatively safe to allow the tourniquet to remain in place until handover to the emergency team at the receiving hospital. It is crucial however to note the time of application. Most purpose-designed tourniquets will have a white tag on which the emergency care provider can note the time of application. Placing a ‘T’ on the forehead of the patient by means of a surgical pen is also a great way of reminding staff that a tourniquet has been applied.

Releasing a tourniquet has its own risks and there are circumstances where removal never makes sense. These later would include limb amputation, shock, and the inability to monitor the wound or continued bleeding. Intermittently releasing them to temporarily restore circulation has been reported to lead to unrecognised, ongoing blood loss and patient death. On a long evacuation, if the conditions seem otherwise safe, waiting 1 hour before attempting release seems like a reasonable time interval. If bleeding starts again, re-secure, note the time and leave it in place.

In conclusion

I believe that tourniquets should be standard issue on all EMS units. Be Safe Paramedical has taken the step of adding a Tactical Tourniquet to all our jump bags as well as our IFAK (Individual First Aid Kits) as a general standard. One simply has to look around at how tourniquets are saving lives to realise that there is no excuse for not having one. This is particular true if one is a professional healthcare provider.

A shark attack which took place in Cape Town a little while ago made headlines around the world because the two elderly gentlemen who rescued the victim, applied a tourniquet to both limbs soon after they got him onto dry sand. There is no doubt that this was the single biggest factor in saving his life. A patient with no blood volume on arrival at a medical facility is inoperable. You cannot replace blood on a scene; however, you can limit its loss and give the patient a fighting chance and the surgeons something to work on.

Give us a call or drop us an email info@be-safe.co.za and you would be extremely surprised at just how cost-effective our Rotation Compression Tourniquet. I challenge you to put one in every jump bag and on every EMS vehicle in your service. I guarantee the results will be phenomenal....

The Responder

“Progress Saves Lives...”

Edition 72 | October 2024

Specials

CritiCare® PneumoPatch™ NV Chest Wound Seal

Product Code: REF:CHEST.WS.PP.NV



R100.00 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

Criticare® Econo Bag Red

Product Code: REF:BAG.ECO.R



R80.50 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

Criticare® RescuPAC EMS Jump Bag

Product Code: REF:BAG.RESCUPAC.CC



R1 250.00 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

CritiCare® PhysioPAC Bag Navy Blue

Product Code: REF:BAG.PHP.CC.N



R230.00 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

The Responder

“Progress Saves Lives...”

Edition 72 | October 2024

CritiCare® AmpPAC Pro Micro Drug Bag

Product Code: REF:BAG.AMP.PRO.M.D



R277.50 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

CritiCare® Rotation Compression Tourniquet

Product Code: REF:CC.R.TOUR



R112.50 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

CritiCare® Grabber Bag (Empty)

Product Code: REF:BAG.GB.CC



R414.20 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

CritiCare® CrashPAC EMS Jump Bag

Product Code: REF:BAG.CP.TB.CC



R2 817.40 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

The Responder

“Progress Saves Lives...”

Edition 72 | October 2024

CritiCare® DrugPAC Medication Bag

Product Code: REF:BAG.DRUGPAC.CC



R1052.25 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

CritiCare® Cardiology Stethoscope

Product Code: STETH.CARDIO



R380.00 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

CritiCare® Rappaport Burgandy Stethoscope

Product Code: STETH.RAPP.BUR



R118.19 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)

CritiCare® CritiBand MkII Trauma Bandage

Product Code: B.Crit.H.100



R122.17 Inclusive of Vat each

Offer valid until 31 October 2024 or while stocks last

[Buy Now](#)